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Table 3.

Getinge is a global provider of innovative solutions for operating rooms, intensive care units, sterilization departments and for life science companies and institutions. Based on our firsthand experience and close partnerships with clinical experts, healthcare professionals and medtech specialists, we are improving the everyday life for people – today and tomorrow.


Flow-i may be pending regulatory approvals to be marketed in your country. Contact your Getinge representative for more information.

The views, opinions and assertions stated by the physician are strictly those of the physician and their practice and do not necessarily reflect the views of Getinge.

Manufacturer - Maquet Critical Care AB · Röntgenvägen 2 SE-171 54 Solna · Sweden · +46 (0)10 335 73 00
US sales contact - Maquet Medical Systems USA · 45 Barbour Pond Drive Wayne, NJ 07470

www.getinge.com

Anesthesia | Ventilation performance

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Printed with permission of the child's parents.

Delivering small volumes with precision

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Case contributed by Dr. Waltraud Bruchelt and Dr. Günter Baumann, Department for Anaesthesiology an Intensive Care Medicine, University Hospital Graz, Austria.

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Table 1. PIP: peak inspiratory pressure

	PIP (mbar)	PEEP (mbar)	MV (l/min)	RR (/min)	TV (ml)	F _I O ₂	pH	p _{cap} O ₂ (mmHg)	p _{cap} CO ₂ (mmHg)	HCO ₃ ⁻ (mmol/l)	BE (mmol/l)	etCO ₂ (mmHg)
Intensive care resp	22	4	0.08	28	2.9	0.35	7.4	30	47	27.1	4.2	n.m.
Flow-i (PCV)	14	5	0.07	24	3	0.31	7.35	48.9	46	26	0.4	34

Table 2.

BE: base excess. cap: capillary. etCO₂: end-tidal CO₂ MV: minute volume. n.m.: not measured RR: respiratory rate. PEEP: positive end expiratory pressure PIP: peak inspiratory pressure PCV: pressure controlled ventilation TV: tidal volume

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