


# Green Tape

## CI 159/.75, CI 159/1, CI 165/.75 & CI 165/1

Green Tape can be used in 121°C gravity sterilizers for 30 minutes or pre-vacuum sterilizers operating at 132°C for 4 minutes and 135°C for 3 minutes. Stripes turn black after exposure to steam sterilization conditions, thus providing identification of processed items.



### Key Features

- 4 minutes at 132°C and 3 minutes at 135°C
- Not made with lead
- Process indicator turns black when exposed to proper steam sterilization conditions

### Performance Characteristics

Indicator stripes change from light green to black

### Storage

The storage conditions are printed on the product packaging

### Shelf Life

36 months from date of manufacture. The expiry date is printed on the product packaging

### FDA 510(K)

K122555

### Standards Classification

- FDA Sterilization Process Indicator
- ANSI/AAMI ST79, Type 1
- ISO 11140-1 Type 1

### SteriTec Part No.

CI 159/.75, CI 165/.75,  
Green Tape 3/4 inch

### Packaging

48 rolls

CI 159/1, CI 165/1,  
Green Tape 1 inch

36 rolls