

# Getinge Clean Enzymatic Detergent

Mild, neutral pH cleaner

Getinge Clean Enzymatic Detergent penetrates and breaks down protein and organic matter.

Getinge Clean Enzymatic detergent combines the power of multi-enzymatic cleaning technology with the advantage of machine spray washing, formulated to remove protein matter from instruments and medical devices.

In addition for usage with automated washers, it is suitable for use in sinks and ultrasonic baths. Enzymatic Detergent is effective in hard and soft water and is suitable for most types of materials.



Available in 5L, 25L and 205L containers

This document is intended to provide information to an international audience outside of the US.

For more detailed information, please contact your local Getinge representative at [getinge.com/contact](http://getinge.com/contact)

[www.getinge.com](http://www.getinge.com)

**GETINGE** 