

Getinge Clean

7900 Sonic Detergent

7900 Sonic Detergent is an effective combination of powerful detergents, wetting agents and emulsifiers.



Application

Excellent for instrument cleaning systems and for removing all types of soils. 7900 Sonic Detergent is ideal for use in ultrasonic cleaning equipment, instrument washers and washer-disinfectors.

Key Features

- Fast, thorough, pollution-free cleaning
- Superior wetting action; thoroughly penetrates all areas, including difficult crevices
- Emulsifies fats and oils, easily removes greasy soils
- Freely rinses with no residue

Usage Information

7900 Sonic Detergent should be used at a dose of 0.2–1.0 oz/gal (2–8 ml/L) depending on water quality. Always follow manufacturer's directions for the particular automatic washer being used. If you have additional questions, consult your Getinge representative.

- Do not add bleach, detergents, or any other chemical to this product. Added chemical may destroy the effectiveness of this cleaner.

NOTE: Some plastics and rubbers are not compatible with detergents or high water temperatures. Check with manufacturer for suitability.

Quality Statement: Confidence in Getinge is the most important quality criteria. This is the hallmark of all our external and internal commitments, activities and products. Products and services supplied by Getinge conform to the agreed terms and expectations. The achievement of these quality goals is the basis for continued competitive and successful enterprise.

Physical and Chemical Properties

- Clear, colorless liquid
- Mild fragrance
- pH value: 10.5–12.8

NOTICE: For institutional and industrial use only. Not for home use. Refer to Safety Data Sheet (SDS) for more information regarding safety, handling and storage and waste disposal methods.

Ordering Information

Getinge Clean	Order Number
4 x 1 Gallon	61301605272
5 Gallons	61301600843
15 Gallons	61301601297

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries • Copyright 2018 Getinge AB or its subsidiaries or affiliates • All rights reserved. • MCV00044134 REVB

Manufacturer: Intercon Chemical Co. 1100 Central Industrial Dr., St. Louis, MO 63110
Sales Office: Getinge USA Sales, LLC, 1777 East Henrietta Road, Rochester, NY 14623

www.getinge.com

