



Thoracic catheter solutions



**DEHP-free
chest tubes**

**Not manufactured
with natural
rubber latex**



PVC thoracic catheters



Soft PVC thoracic catheters

- Advanced taper tip designed for easy placement
- Soft, thermo-conforming PVC
- 2 cm depth markings aid in accurate catheter placement
- Gentle spiral hole configuration
- Tapered oval drainage eyes
- Blue radiopaque stripe confirms placement upon X-ray

Firm PVC thoracic catheters

- Open tip designed for fast, secure connection to chest drain tubing
- Smooth, atraumatic tip



Blue radiopaque line confirms placement upon x-ray



Silicone thoracic catheters

Round silicone catheters

- Soft silicone designed to conform to anatomy for improved patient comfort
- Silicone provides drainage visibility
- Blue radiopaque stripe confirms catheter placement upon X-ray

Mediastinal catheters

- Low profile design
- Choice of eyelet configurations

Extended length (XL) catheters

- Extended length catheter with 15 to 18 eyelets
- Low profile 20Fr and 24Fr sizes

Blue radiopaque line confirms placement upon x-ray



Indication by design

Getinge offers a full line of catheter products designed to meet the needs of cardiac, thoracic, and pediatric patients. All of our thoracic catheters are DEHP-free and not made with natural rubber latex.

Getinge also offers a variety of accessories to secure any size thoracic catheter to adult, pediatric or mobile chest drains.

Your application. Your choice.

Soft PVC thoracic catheters

Straight catheter

Size (Fr)	Description	Length (in)	# Eyelets	Code #
8	Straight	8	2	8008
10	Straight	8	2	8010
12	Straight	16	4	8012
16	Straight	16	4	8016
20	Straight	22	4	8020
24	Straight	22	5	8024
28	Straight	23	6	8028
32	Straight	23	6	8032
36	Straight	23	6	8036
40	Straight	23	6	8040

Right angle catheter

Size (Fr)	Description	Length (in)	# Eyelets	Code #
12	Right angle	16	4	8112
16	Right angle	16	4	8116
20	Right angle	22	4	8120
24	Right angle	22	4	8124
28	Right angle	23	5	8128
32	Right angle	23	5	8132
36	Right angle	23	5	8136
40	Right angle	23	5	8140

Firm PVC thoracic catheters

Straight catheter

Size (Fr)	Description	Length (in)	# Eyelets	Code #
12	Straight	14	4	15012
16	Straight	14	4	15016
20	Straight	19	4	15020
24	Straight	19	5	15024
28	Straight	20	6	15028
32	Straight	20	6	15032
36	Straight	20	6	15036
40	Straight	20	6	15040

10 units per case. Please consult your local representative for availability.

Silicone thoracic catheters

Straight catheter

Size (Fr)	Description	Length (in)	# Eyelets	Code #
16	Straight	20	4	14016
20	Straight	20	5	14020
24	Straight	20	6	14024
28	Straight	20	6	14028
32	Straight	20	6	14032
36	Straight	20	6	14036

Extended length (XL) catheter

Size (Fr)	Description	Length (in)	# Eyelets	Code #
20	Extended length straight	26	15	14820
24	Extended length straight	26	18	14824

Right angle catheter

Size (Fr)	Description	Length (in)	# Eyelets	Code #
16	Right angle	20	4	14216
20	Right angle	21	5	14220
24	Right angle	21	6	14224
28	Right angle	21	6	14228
32	Right angle	21	6	14232
36	Right angle	21	6	14236

Mediastinal catheter (Oval)

Size (Fr)	Description	Length (in)	# Eyelets	Code #
7	Mediastinal, eyes - top & bottom	20	6	14407
9	Mediastinal, eyes - top & bottom	20	9	14409
11	Mediastinal, eyes - top & bottom	20	9	14411
7	Mediastinal, eyes on top	20	4	14607
9	Mediastinal, eyes on top	20	6	14609
7	Mediastinal, eyes on sides	20	5	14807
9	Mediastinal, eyes on sides	20	5	14809

Right angle catheter

Size (Fr)	Description	Length (in)	# Eyelets	Code #
24	Right angle	19	4	15224
28	Right angle	20	5	15228
32	Right angle	20	5	15232
36	Right angle	20	5	15236



Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries. · Copyright 2020 Getinge AB or its subsidiaries or affiliates · All rights reserved · ⚠ CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician. Refer to Instructions for Use for current indications, warnings, contraindications, and precautions.

Sales Office · Getinge USA Sales LLC · 1 Geoffrey Way · Wayne, NJ 07470

www.getinge.com

MCV00046186 REVE
EXP08/22