

Wash-Checks H Holder

WC 106

The Wash-Checks H Holder is an auxiliary device for routine monitoring of automated wash processes for hollow instruments in washer-disinfectors and ultrasonic washers.

Wash-Checks H Holder should only be used with Wash-Checks H indicators (in washer-disinfectors) or Wash-Checks UH indicators (in ultrasonic washers). Wash-Checks H or Wash-Checks UH is inserted into the Wash-Checks H Holder, simulating contaminated rigid lumen instruments.

Key Features

- Reusable holder
- Electro polished stainless steel AISI 304
- Use with WC 105 and WC 110

Performance Characteristics

The Wash-Checks H or Wash-Checks UH is inserted into the Wash-Checks H Holder. The holder is then attached to a luer lock connector located on the washer rack.



Physical Characteristics

- **Dimensions:** 197 mm X 11 mm
- **Weight:** 55 g

Shelf Life

Not applicable

SteriTec Part No.

WC 106, Wash-Checks H Holder

Packaging

Each