

Wash-Checks Pro Protein Detection Kit 6" Swab WC 115

Wash-Checks Pro is a quick and easy test used to detect residual protein in hollow lumen surgical instruments.



Key Features

- 6 inch swab (152 mm)
- The higher level of contamination of protein the darker and faster the blue result
- Accuracy of 1 µg in less than 10 seconds
- No incubation required

Performance Characteristics

This protein test can detect protein residues within 1 µg sensitivity and 10 seconds. When protein is present, the solution will turn from brown to blue. This test provides a semi-quantitative measure of surface cleanliness; the higher the level of contamination of protein (bioburden), the darker and faster the color changes to blue.

Storage

Storage conditions are printed on the product packaging. Keep away from light and heat.

Shelf Life

2 years from manufacture of reagent when refrigerated and 6 months when not refrigerated.

Standards Classification

Consistent with ISO 15883, routine monitoring of washer disinfectors

SteriTec Part No.

WC 115, Wash-Checks Pro

Packaging

2 kits of 25 tests & 25 swabs

FTP-350-1 Rev. B



SteriTec Products Mfg. Co., Inc.,
Member of the Getinge Group
74 Inverness Drive East
Englewood, CO 80112 USA

www.getinge.com/steritec