

Getinge Clean

Renuzyme Foam Spray

Renuzyme Foam Spray is a ready-to-use dual enzyme pre-treatment foam designed for point-of-use, reprocessing delays, and in the Sterile Processing Department to enhance the manual cleaning process.



Application

For point-of-use (prior to transport) to begin decontamination process as soon as possible, post procedure per AAMI ST79 and AORN recommendations.

Key Features

- Improves effectiveness of cleaning process by increasing product contact time with soiled surface
- Prevents drying of blood and other bioburden on surfaces before decontamination
- Effective at room temperature
- Will not stain or damage glass, metal, plastic or equipment
- Pleasant, green apple fragrance

Usage Information

Please consult Instructions for Use (IFU) for specific usage information. Device manufacturer IFU should be consulted for cleaning instructions.

- Do not add bleach, detergents, or any other chemical to this product. Added chemical may destroy the effectiveness of this cleaner.

Quality Statement: Confidence in Getinge is the most important quality criteria. This is the hallmark of all our external and internal commitments, activities and products. Products and services supplied by Getinge conform to the agreed terms and expectations. The achievement of these quality goals is the basis for continued competitive and successful enterprise.

Physical and Chemical Properties

- Clear, light amber color
- Green apple fragrance
- pH value: 9-10

NOTICE: For institutional and industrial use only. Not for home use. Refer to Safety Data Sheet (SDS) for more information regarding safety, handling, storage and waste disposal methods.

Ordering Information

Getinge Clean

Order Number

12 1-Quart Bottles

61301604584

5 Gallons

61301606712

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries • Copyright 2017 Getinge AB or its subsidiaries or affiliates • All rights reserved. • MCV00044126 REVA

Manufacturer: Intercon Chemical Co. 1100 Central Industrial Dr., St. Louis, MO 63110

Sales Office: Getinge, 1777 East Henrietta Road, Rochester, NY 14623

www.getinge.com

GETINGE 