

# Getinge Clean

## Stainless Steel Chamber Cleaner

Stainless Steel Chamber Cleaner is an acid detergent for use in cleaning and brightening chambers and exteriors of equipment such as sterilizers and washer-disinfectors.



### Application

For use in removing scale, discoloration and mineral deposits from chambers and washer-disinfectors including: cart, cage and rack washers.

### Key Features

- Strong dissolving action on scale and discoloration from mineral deposits such as calcium and magnesium carbonate
- Free rinsing
- Complete solubility in water

### Usage Information

This product may be used as-is, or diluted up to 1:1 with water.

- Do not add bleach, detergents, or any other chemical to this product. Added chemical may destroy the effectiveness of this cleaner.

Please consult Instructions for Use (IFU) for dilution rates and specific information. Equipment manufacturer IFU should be consulted for cleaning instructions.

**Quality Statement:** Confidence in Getinge is the most important quality criteria. This is the hallmark of all our external and internal commitments, activities and products. Products and services supplied by Getinge conform to the agreed terms and expectations. The achievement of these quality goals is the basis for continued competitive and successful enterprise.

## Physical and Chemical Properties

- Clear, colorless thin liquid
- Characteristic fragrance: no fragrance added
- No added dyes
- pH value: <1.0-1.5

Notice: For institutional and industrial use only. Not for home use. Refer to SDS for more information regarding safety, handling and storage and waste disposal methods.

---

## Ordering Information

### Getinge Clean

### Order Number

4 x 1 Gallon

61301605284

---

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries • Copyright 2017 Getinge AB or its subsidiaries or affiliates • All rights reserved. • MCV00044181 REVA

**Manufacturer:** Intercon Chemical Co, 1100 Central Industrial Dr., St. Louis, MO 63110  
**Sales Office:** Getinge, 1777 East Henrietta Road, Rochester, NY 14623

[www.getinge.com](http://www.getinge.com)

**GETINGE** 