



SERVINO – iNO delivery system

Delivering therapy when it is critical



Ready

Preparing

On-screen help, with images and text, makes for a quick and easy setup, allowing for a faster response to your patient.

Go

Running

You can deliver a clinically relevant NO dose, and NO, NO₂ and O₂ levels are clearly displayed and trended for 72 hours for careful monitoring.

Free

Weaning

The controlled, monitored reduction of NO levels helps you safely wean your patient, reducing the risk of the rebound phenomenon and harm.

Delivering therapy when it is critical

While responding to a patient in distress, you need to react quickly. The Servino™ device is built to provide precision and accuracy to your NO delivery, enabling you to focus on your patient.

Delivering safe treatment is vital, that's why we have a highly advanced NO delivery and monitoring device. It delivers nitric oxide (NO) – a vasodilator – in harmony with many validated ventilators to improve oxygenation, trying to prevent more invasive treatment.

Servino is safeguarded to supply a clinically relevant dose. The sensitive sensors allow the NO gas to have an accurate and continuous delivery, reducing the risk of rebound pulmonary hypertension and patient harm.

Servino also helps decrease your workload with its easy-to-use interface, quick startup for rapid response to babies in crises, automatic gas sensor adjustment and a simple-to-use gas supply calculator.





Getinge is a global provider of innovative solutions for operating rooms, intensive care units, sterilization departments and for life science companies and institutions. Based on our firsthand experience and close partnerships with clinical experts, healthcare professionals and medtech specialists, we are improving the everyday life for people – today and tomorrow.

This document is intended to provide information to an international audience outside of the US. Servino may be pending regulatory approvals to be marketed in your country. Contact your Getinge representative for more information.

Manufacturer · Maquet Critical Care AB · Röntgenvägen 2 SE-171 54 Solna · Sweden · +46 (0)10 335 73 00

www.getinge.com