

Cardiothoracic Surgery:

A full complement of advanced chest drainage products



Atrium Ocean series wet suction, water seal chest drains use a familiar water seal operating system with traditional wet suction for all thoracic drainage applications. Available in single, dual, ATS and infant.



Atrium Oasis series dry suction, water seal chest drains use a familiar water seal operating system with dry suction regulation for all thoracic drainage applications. Available in single, dual, ATS and infant.



Atrium Express series dry suction, dry seal chest drains use a dry seal operating system with dry suction regulation for all thoracic drainage applications. Available in single and ATS.



Getinge's mobile drainage technology includes the dry suction, dry seal **Express Mini 500 Chest Drain** for air and fluid drainage as well as the **Pneumostat Chest Drain Valve** for air-only management.



Getinge's **PleuraGuide Disposable Chest Tube Kit** has a comprehensive assortment of everything likely needed for a bedside chest tube insertion.

Cardiothoracic Surgery: A full complement of advanced chest drainage products



PVC Catheters: Getinge's full line of PVC catheters include soft (8000) and firm (15000) PVC in straight and right angle configurations with sizes from 8Fr to 40Fr to accommodate a wide variety of thoracic drainage applications.

DEHP
Free



Silicone Catheters: Getinge's full line of silicone thoracic catheters includes round, extended length (XL) and flat mediastinal configurations to accommodate a wide variety of thoracic drainage applications.



Blood Bags: Autotransfusion blood bags are used in conjunction with Getinge Chest Drains to provide reliable and cost-effective blood management solutions.



Accessories: Include Y-connectors, step-down connectors, tubing packs and more to attach to any Getinge Chest Drain.

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries • Copyright 2017 Getinge AB or its subsidiaries or affiliates • All rights reserved • ⚠ CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician. Refer to Instructions for Use for current indications, warnings, contraindications, and precautions.

www.getinge.com

GETINGE 

MCV00047427 REV B
EXP07/19